Check Stab Instruments

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Functions of the Check Stab α 2015 iDry

- *Quick and easy measurements of tartaric acid stability and calcium saturation using electrochemical conductivity measurements at various conditions.*
- *Mini Contact: Determines KHT precipitation expressed in μ Siemens and %. This value indicates the relative stability of the wine sample being tested.*
- *Mini contact with forecast : Determines KHT precipitation expressed in μ Siemens and % It uses CheckFuture® to predict future knock down of crystals.*
- *Saturation Point: Automatically, the analysis calculates the saturation temperature of the test sample expressed in °C. It has also the possibility to evaluate the data with the total polyphenols index.*
- *Saturation Point of Calcium: This analysis calculates automatically the saturation temperature value of calcium.*
- *Special Functions: ISS®, ISTC50.*

Check Stab α 2015 iDry uses a more evolved and innovative technology for the control of tartaric stability. **Check Stab α 2015 iDry** has a new refrigeration system for the sample without the use of liquid glycol; this system allows the instrument to shorten the duration of the analysis, be more precise within the measuring system and at the same time be more ecological. The instrument is equipped with a 9,9" touch screen colour display, interfaced with the software CheckStab.NET® by means of LAN, meaning it can be connected directly to the winery's computer network. **Check Stab α 2015 iDry** is a completely automatic analyzer and can satisfy the needs of all wineries, from small companies to big industry, thanks to the possibility to modulate a system of up to 16 units, all managed by one software, Check Stab NET®, storing data on the same Dbase and an automatic KHT dispenser.

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View of the automatic KHT dispenser

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Some characteristics of the Check Stab α 2015 iDry

Possibility to insert data of wine sample and start the analysis both at the instrument or from the PC with synchronization of the dbase. The instrument stores the last 10 analyses. Possibility to connect to Bar Code reader to insert sample wine data. Possibility to connect more than one instrument together (max 16 units), with each unit doing an independent analysis. Automatic dispensing system for adding KHT without the intervention of the lab technician.

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
*Acquiring and recording measurement data is becoming more and more important. **CheckStab Instruments®** has developed a new analytical software system to manage data. This permits the lab to study, store, and comment on an analysis with extreme simplicity. The entire system is managed by use of a personal computer.*

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What Check Stab α2015 iDry and Check Stab.NET can do:

- Visualize measurement data (conductivity, wine temperature, bath temperature, pre-set temperature, etc.).
- Visualize graphics in real time
- Possibility to start the analysis directly from instrument, and automatic synchronization with the PC.
- Visualize user status, using incorporated graphic symbols (state of cooler; heater; probe positions, etc.).
- Memorize the analysis on the data-base, with the possibility to export to a text file.
- Possibility to get from Check Stab analyses done, to visualize results in graphic form, to visualize results in text form, to save on a disk, export reports in various formats (PDF, Excel, Text) or print the graphics using a printer.
- Visualize the saturation analysis in a normal graphics mode or in a differential mode.
- Visualize the mini contact analysis in absolute mode or in relative mode.
- Visualize the analysis in single mode or multiple mode, up to a maximum of 10 analysis contemporaneously.
- Interface for LIMS, sends automatically reports by e-mail.

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Data given in this brochure are indicative and can vary at the discretion of the manufacturer at any time and without obligation or notice.

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